

THE NEW EVO



TAKE ANOTHER LOOK...

For 40 years, we've manufactured animal-specific ultrasound right here in the US. Now we've done it again... the next generation of our flagship system, the **NEW EVO**.

With new high-resolution probes, beautiful images and improved battery performance, the **EVO III** provides high-quality ultrasound in a rugged, portable machine.

- New broadband probes for all species, including the all-new hybrid **eCli6**
- Completely redesigned image software with the latest technology
- Customized, veterinary-specific exam presets
- Sunlight-readable LED display
- WiFi streaming via EVOStream app
- Ruggedized case, splash resistant, and portable [just over 6 lbs]—sealed keyboard and probe connections allow easy cleaning
- Interchangeable batteries, for continuous use
- Designed and manufactured in Colorado

E.I. MEDICAL IMAGING
1.866.365.6596 • www.eimedical.com/the_new_evo



BROADBAND TRANSDUCERS



eCL3
2-5 MHz macro convex



eCLi4
3-7 MHz curved rectal



eCLi6
3-9.5 MHz hybrid curved linear **NEW!**



eC6
4-9 MHz micro convex



eC9
5-10 MHz micro convex



eL7
5-9 MHz linear rectal



eL14
6-14 MHz tendon



L3ASE
Linear transducer,
Bovine carcass merit.



eC90PU
5-10 MHz broadband array for transvaginal oocyte aspiration in horses and cattle. Adjustable handle, small height, lightweight.



PROBE ACCESSORY



I.C.E.®
The IBEX Customizable Extension lessens fatigue and stress on the user's arm and shoulder.
New, improved version is 25% lighter, smaller diameter, more durable than previous models.
The I.C.E. extension allows easy insertion or removal of the probe without damaging cables.
Compatible with all IBEX and EVO rectal probes.

A WORD ABOUT OUR CABLES...

IBEX® cables tested to 1.2 million flex cycles without breaking

- Custom designed cables are abrasion resistant
- Only manufacturer to design all transducers specifically for animal use
- Kevlar reinforced

...AND OUR NEW TRANSDUCERS

As part of our commitment to innovation, we have re-designed our probes for improved penetration and resolution. Broadband, variable frequency, our transducers are designed specifically for rugged environments. But if the worst happens, they can be repaired right here in our Colorado facility.

E.I. MEDICAL IMAGING

E.I. Medical Imaging® was founded in 1984 with the vision of developing state-of-the-art, real-time ultrasound to serve veterinarians and livestock producers world-wide.

Our products evolved with the needs of the market. We are the only manufacturer of portable ultrasound systems engineered and serviced in the United States specifically for the animal industry.



CONTACT US

1.866.365.6596 [toll free]
1.970.669.1793
1.970.669.1902 [fax]
info@eimedical.com
www.eimedical.com

